

Εύρηκα

研究 CTEN 在大腸 癌細胞中的影響

報告人: 江祐辰
指導老師: 魏宏仁
指導教授: 廖憶純



Εύρηκα

背景



研究動機

Eύρηκα


- 台灣是個美食天堂
- 外食頻繁
- 大腸癌是台灣罹患率最高的癌症之一




CT

- C 集中胞附要角


- CT 過程



(a)



(c)



(b)

K α

質，
於細
寅重
生。

遷移的

β

文獻

Eύρηκα

- 正常情況下 CTEN 表現量低(集中附著點)
- 大量 (表現量高)的 CTEN 也在大腸癌細胞的細胞核中被發現，這是正常細胞中沒有被觀察到的現象
- CTEN 在大腸癌細胞核中可與 beta-catenin 結合並促進腫瘤發生。



Εύρηκα

目標



目標

Εύρηκα

- 瞭解位於細胞核的 CTEN 在腫瘤發生過程中所扮演的角色，進而闡明 CTEN 促進腫瘤發生的分子機轉，預期可進而提供發展癌症治療的新標的。



Εύρηκα

方法




染色質免疫沉澱

Εύρηκα


- Chromatin Immunoprecipitation
- 簡稱ChIP
- 功能：確定特定蛋白是否結合特定基因組區域



Εύρηκα



Εύρηκα



步驟(Crosslink)

Εύρηκα

1

- PBS潤洗細胞
- 加入 DMEM & 甲醛 開始 Crosslink

2

- 加入 甘胺酸 停止 Crosslink
- PBS潤洗細胞

3

- 加入PBS將細胞刮入離心管，離心後丟掉上清液[+C]
- 加入PBS，分裝到微量離心管+C

步驟(Shearing)

Εύρηκα

1

- 以 Buffer A & NP40 [P.I.] 懸浮細胞(震盪)+C
- 以 MNase rxn buffer 懸浮。加入 Mnase，37°C乾浴

2

- 加入 MNase stop buffer+C
- 以 RIPA Lysis Buffer 懸浮

3

- 超音波剪切
- 離心後收集上清液* (取約1/4為Input)

步驟(wash&block beads)

1

- 震盪agarose beads，分裝至微量離心管
- 加dilution buffer，離心後置冰上一分鐘，丟掉上清液×2

2

- 加入BSA(小牛血清白蛋白)&dilution buffer
- Rotate for 1hr (4°C)

3

- 加dilution buffer，離心後置冰上一分鐘，丟掉上清液
- 加dilution buffer 沖掉多餘的bead& 重新懸浮

步驟(Pre-clearing)

Εύρηκα

1

- 加 dilution buffer in [P.I.]& BSA-blocked beads
- 4°C水浴 1 hr 後，離心，**收集上清液**到新的微量離心管

2

- 加入抗體，4°C水浴 16 hr
- 加入 blocked beads，4°C rotate 2 hr

3

- 離心後置冰上一分鐘
- 以 low-salt wash buffer 沖洗(緩慢pipet)+C

步驟(Pre-clearing)

Εύρηκα

4

- 將混合液裝至新的微量離心管+C
- 以high salt wash buffer 沖洗+C

5


- 以 LiCl wash buffer 沖洗+C
- 以 TE Buffer 沖洗，裝至新微量離心管+C

6

- 加 ChIP elution buffer，65°C水浴搖晃半小時
- 離心，將上清液裝到新的微量離心管

步驟(De-crosslinking)

- 
- 融化Input，加入elution buffer
 - 加 5M NaCl到每個 sample，65°C水浴反應

- 
- 加 Tris-HCl(ph6.8) & EDTA
 - 加 Protein K

- 
- 65°C水浴
 - 純化DNA

步驟 (Crosslink)

Εύρηκα

- PBS rinse cell
- Refill 9ml DMEM+270 μ l 37% formaldehyde, shake 10min at regular temperature(@RT)
- Add 1ml 1.25M Glycine, shake 5min @RT
- Wash with icy PBS \times 2
- Add 3ml icy PBS to each plate and scrape it into 15ml falcon tube
- Centrifuge (1000g, 5min @4°C), discard supernatant
- Suspended with PBS, transfer 1ml to each eppendorf



步驟 (Shearing)

Εύρηκα

- Centrifuge (2000g, 5min @4°C), discard supernatant
- Suspended with [100μl BufferA+7μl **NP40**] ([P.I.]
- Vortex, rest on ice for 2min, vortex
- Centrifuge (2000g, 5min @4°C), discard supernatant
- Suspended with 360μl MNase rxn buffer
- Add 1ml **MNase**, dry bath @37°C, flip every 5min
- Add 80μl MNase stop buffer



步驟 (Shearing)

Εύρηκα

- Centrifuge (9000g, 5min @4°C), discard supernatant
- Suspended with 400μl **RIPA lysis Buffer**×[P.I.]
- Sonicate:[10W, on 30sec ; off 30sec] ×4
- Centrifuge (8000rpm, 10min @4°C)
- **Collect supernatant***, take aside approximately 1/4 as "Input"
- Store @80°C or move on



步驟 (Wash & Block Beads)

(do it 2hr beforehand)

- Vortex agarose beads until fully suspended
- Transfer 25 μ l beads to each eppendorf, dilute to 1ml by dilution buffer
- Centrifuge (2000g, 1min @4°C), let beads settle on ice for 1min, discard supernatant
- Add 1ml dilution buffer, invert tubes 5 times
- Centrifuge (2000g, 1min @4°C), let beads settle on ice for 1min, discard supernatant



步驟 (Wash&Block Beads)

- Add 100 μ l 10mg/ml **BSA** & 900 μ l dilution buffer
- Rotate for 1hr @4°C
- Centrifuge (2000g, 1min @4°C), let beads settle on ice for 1min, discard supernatant
- Wash beads once with 1ml dilution buffer
- Suspended with 25 μ l dilution buffer



步驟 (Pre-clearing)

Εύρηκα

- Add 900μl dilution buffer(P.I.)
- Add 25μl BSA-blocked beads
- Incubate @4°C F-wheel for 1hr
- Centrifuge (7500rpm, 3min @4°C), collect supernatant to other eppendorf
- Add **antibodies** to each eppendorf, incubate 16hr @4°C
- Add 50μl blocked beads, rotate for 2hr @4°C
- Centrifuge (2000g, 1min @4°C), let beads settle on ice for 1min



步驟

Εύρηκα

- Wash with high salt wash buffer, centrifuge, discard supernatant
- Wash with 1ml LiCl wash buffer, centrifuge, discard supernatant
- Wash with 1ml TE wash buffer, transfer mixture to the eppendorf, centrifuge, discard supernatant
- Add 210µl CHIP elution buffer, incubate shaking for 30min @65°C, centrifuge(11000rpm, 3min @4°C)
- Transfer 200µl of supernatant to the new eppendorf



步驟 (De-crosslinking)

Εύρηκα


- Thaw input and add elution buffer to 200μl
- Add 8μl of 5M NaCl to each sample, incubate, shake at 65°C O/N
- Add 20μl 1M Tris-HCl (pH6.8) & 20μl 0.25M EDTA
- Add 2μl Proteinase K(20mg/ml)
- Incubate and shake for 1hr @60°C
- Purify DNA (kit: GP1002)



跑

• 以 D

	DHX32-1	DHX32-2	DHX32-3	Average	Average*	Delta-Ct	fold enrichment
Input (5%)	21.67	21.51	20.91	21.36333333	17.04140524		
RNAP	28.24	28.28	26.1	27.54	27.54	2.256666667	0.209254903
IgG	25.43	25.27	25.15	25.28333333	25.28333333	0	1
CTEN	23.64	25.11	24.29	24.34666667	24.34666667	-0.936666667	1.914100614
Input (5%)							
RNAP						-1.02	2.02791896
IgG						0	1
CTEN						1.176666667	9.042155114
Input (100%)							
RNAP						-0.39	1.310393404
IgG						0	1
CTEN						-1.83	3.555370725
STDEV	RNAP						
	0.5						
Average	RNAP					2.02791896	1.310393404
	1.1					1	1
TTEST	RNAP vs IgG		CTEN vs IgG			1.914100614	9.042155114
	0.747395234		0.149606743				3.555370725



K α

B

目前結論

Εύρηκα

跑基因定序後發現

Cten 會和與細胞週期、移動有關的基因接合。



未來展望

Εύρηκα

- 釐清運送 CTEN 至細胞核的機制
- 細胞核 CTEN 與促進腫瘤發生的關連性。



參考資料

Εύρηκα

- <https://www.ncbi.nlm.nih.gov/pubmed/19214987>
- https://en.wikipedia.org/wiki/Chromatin_immunoprecipitation
- <https://ndltd.ncl.edu.tw/cgi-bin/g32/gswweb.cgi/login?o=dnclcdr&s=id=%22102NTU05106070%22.&searchmode=basic>
- <https://ndltd.ncl.edu.tw/cgi-bin/g32/gswweb.cgi/login?o=dnclcdr&s=id=%22102NTU05106070%22.&searchmode=basic>
- <https://www.thermofisher.com/order/catalog/product/88216>
- <https://www.viogene.com/products/39>



Εύρηκα

最後感謝



Εύρηκα

- 廖憶純教授
- 實驗室的所有學長姐
- 魏宏仁老師
- 生專的同學
- 所有聆聽的人



Εύρηκα

Thank You

